

MEYRA



MC2-RS

+ THE POWERFUL MC2
GOES NEXT LEVEL IN STYLE



Crash tested for safety
Complies with international
ISO safety standards for
crash safety

+ **Stylish, sporty and powerful**
Awesome red and black colour scheme
and Meyra RS or Recaro seating

+ **High quality German design
and manufacture**
Made to last with the highest quality
components

+ **Heavy duty**
User weight up to 160kgs

+ **Highly manoeuvrable**
Conquer the indoors and
outdoors with ease

+ **Crash tested**
User can be transported by
bus and maxi taxi

Grunt and good looks

RIDE IN STYLE


MC2-RS: Power and good looks, together at last. This special edition of the beautifully engineered MC2 takes you to the next level in style and power, with an MC2-RS exclusive seating system, 15" black anodised drive rims, optional wide tyres for outdoors and the grunt to match, as well as customisation and accessories galore!




15" drive wheels



Side guard lighting with custom engraving



Black-red RS design



Recaro seat in black leatherette

TECHNICAL DATA

MC2-RS

| | | | |
|-----------------------------------|--|------------------------|---|
| Seat width | 38 - 65 cms | Turning radius | Approx. 84 cms |
| Seat depth | 40 - 56 cms | Driving range | 35 - 40 kms with 73 ah batteries ¹ |
| Backrest height | 45 - 50 cms* *depending on system selected | Chair weight | Approx. 100 kgs, no batteries |
| Legrest length | 35 - 50 cms | Maximum client weight | 160 kgs |
| Armrest height | 18 - 29 cms | Motor power rating | 350 watt 4 pole German motors |
| Seat height - front | 53 - 65 cms *depending on tyre & seating combinations | Maximum gradient | 15 degrees |
| Overall length - with legrests | Approx. 113 cms | Mechanical seat tilt | 0 to 10 degrees |
| Overall length - without legrests | Approx. 83 cms | Power seat tilt | 0 to 50 degrees, tilt in space |
| Overall width | 63.5 cms | Seat raise | 300mm |
| Wheel size front | 260 x 70mm | Mechanical backrest | -10 to + 30 degrees |
| Wheel size rear | 380 x 75mm | Power backrest recline | -10 to + 50 degrees |
| Maximum obstacle climbing height | Approx. 60mm | Compliance testing | ISO7176 strength & durability standards ISO7176-19 crash test requirements |

¹ Distance will be affected by client weight, terrain and battery condition

Dimension tolerances
± 10 mm ± 2°

Data subject to changes

MEYRA
AUSTRALIA



MEYRA Australia
5 Boulder Rd Malaga
Western Australia 6090



Tel (08) 9248 3560
Fax (08) 9209 1999



sales@meyraaustralia.com.au
www.meyraaustralia.com.au