

MEYRA

XL & XXL

+ TOUGH
AND COMFORTABLE



User weight up to 250 kg

*The perfect bariatric wheelchair
—up to 200 kg for XL and up to 250 kg for
XXL*



Rear wheel drive

*The classic stable drive configuration
—great for indoors and outdoors*



Built tough with lots of options

*Reinforced components for strength
and stability, with a wide range
of standard and power options*



LED lighting system

*Safely light your way
with long-lasting LEDs*

The chair for rugged types

TOUGH AND STABLE


The Meyra XL and XXL don't back down. Built to move users up to 200 kg (XL) or 250 kg (XXL), with a great stable and sturdy design, the XL and XXL are tough. Thanks to the power tilt-in-space and recline function, as well as the choice of powered separate or centre mounted footrests, you can be tough and comfortable.



Reinforced and built tough



Powerful 2 x 300W or
2 x 350W motors



Powered tilt-in-space



Long-lasting
LED lighting system

TECHNICAL DATA

XL & XXL

Seat width	530 - 700
Seat depth	430 - 600
Back rest height	530 - 570
Leg length adjustment	370 - 510
Armrest height adjustment	170 - 280
Front seat height options (with /without cushion)	420 - 470 / 460 - 560 / 530 - 630
Length with footrests	1,150
Length without footrests	850
Width overall	670
Height without headrest	1,020 - 1,150
Transport length	850
Transport height minimum	650
Front wheel size	260 x 70 (10" x 2 ¾")
Rear wheel size	356 x 75 (14" x 3")

Max. obstacle height	60
Turning radius	950
Driving range with 70ah batteries	30 - 35
Empty transport weight min.	approx. 130 kg
Client weight max XL / XXL	200 kg / 250 kg
Driving speeds	6 kmh / 10 kmh
Motor options	2 x 300W / 2 x 350W
Driving range max.	35 km
Permitted gradient	13 %
Electric tilt-in-space	0° to +18°
Mechanical tilt-in-space	0° to +10°
Power backrest recline	-10° to +50°
Mechanical backrest recline	-10° to +30°

Dimensions in mm, Weights in kg, Speed in km/h, Tolerance ± 10 mm